



Restylane
GALDERMA



آفامید
آفتاب فروزان ایرانیان
مدیریت دارو

تهران، خیابان مطهری، میرزای شیرازی شمالی
خیابان ناصید، پلاک ۲۲، طبقه دوم تلفن: ۰۲۱ ۸۸ ۱۰۲ ۱۰۲

www.afamed.org

[afamed.co](https://www.instagram.com/afamed.co)

Product Catalogue
2020 - 2021

Restylane is the most studied wrinkle filler in history and has had dozens of articles about its development, uses and advancements have been published in top medical journals.

Development of Restylane

Restylane was studied and researched for many years before being introduced to the public for use in filling as well as augmentation. It was first used in Europe in 1996 and soon crossed over to the United States, as well as to about 80 other countries. Not only is Restylane the number one most studied wrinkle filler, but it is also the number one used wrinkle filler and augmentation injectable gel.

In 2004, Restylane began being used in the US, and millions of procedures have been done using the clear injectable gel for the purpose of smoothing wrinkles and also in augmentation.

Restylane is the world's broadest range of hyaluronic acid (HA) dermal fillers. Based on two innovative technologies (NASHA & OBT)* that restore your facial appearance.

Restylane dermal fillers offer a versatile solution to address your individual needs, in a number of facial areas. The clear gel, made from hyaluronic acid, has been FDA approved for wrinkle filling and then it was approved to be used in the deep dermis. This opens the doors for even more advancements in cosmetic procedures and longer lasting results.

With over 40 million treatments worldwide and more than 24 years of clinical experience, Restylane products are recognized and trusted to provide documented safe, instantly visible and natural-looking results that can last up to 12-18 months, depending on the treatment area.

***NASHA:** Non Animal Stabilized Hyaluronic Acid (HA particles were clearly emphasized and ballon shaped rather than a round pearl. The process uses bacterial fermentation to obtain a high degree of purity of natural HA. Potentially harmful components, such as proteins and endotoxins are minimised, and those of animal origin in particular, are excluded by the process. The natural HA molecule forms a highly viscous gel due to the natural entanglement of the long chains of the molecule.)

***OBT:** Optimal Balance Technology(This technology offers a broad range of gel textures, all containing the same concentration of Hyaluronic Acid gels (HA), which means that the products can serve a wide range of functions. Softer gels distribute more evenly to smooth lines and wrinkles, and firmer gels can deliver more lift and volume.)

Injectable Gel With Lidocaine

DEFYNE™

The EMERVEL Collection

What is Restylane Defyne with Lidocaine ?

Restylane Defyne is a dermal filler injected under the facial skin. It is made out of hyaluronic acid – a naturally occurring sugar, found in the body that gives the skin moisture, volume, and elasticity – and a small amount of lidocaine, which is a commonly used local anesthetic to numb the skin.

Restylane Defyne is made using a manufacturing process called Optimal Balance Technology (OBT), which creates a smooth, injectable gel that can give your skin a natural, soft look.

Restylane Defyne is the only filler on the market that has been shown in a clinical study to help maintain movement all expression.

CLINICAL STUDIES:

Restylane Defyne was tested to ensure that it worked properly (effective) and was safe to use, in 162 patients before it was approved for use in all patients. These patients included 156 women, 6 men, and were from various ethnic groups with multiple skin tones (pale to dark skin). Each of these patients expressed desire to improve the appearance of their facial wrinkles and folds (laugh lines). The study was created to compare two different dermal filler gel products. It was conducted by 11 different doctors in 11 different locations across the United States. Patients who participated in the study received an initial treatment with Restylane Defyne Touch-up injections were given as needed, and each patient was monitored for 12 months following the procedure to measure the effects of the treatments.

For each patient, the amount of product used was based on their individual needs. In order to objectively measure improvements in appearance over 18 months, each patient was evaluated by a doctor other than one who had performed the treatment. Patients also evaluated themselves. To measure side effects following the procedure, each patient kept a diary for 21 days to record side effects. These were then shared with the doctors.

The study doctors reported:

93% of subjects had at least a 1-point improvement in wrinkle appearance 6 months after treatment.
87% of subjects had at least a 1-point improvement in wrinkle appearance 1 year after treatment.



Defyne
your
Expression





VOLYME
Can Restore Your
Youthful Appearance

Injectable Gel With Lidocaine

VOLYME[™]

The EMERVEL Collection

What is Restylane volyme with Lidocaine?

Restylane Volyme Lidocaine is a Hyaluronic Acid based volumizing dermal filler designed to create a volume effect when injected into the supraperiosteal zone or subcutaneous tissue. The filler is based on the patented Optimal Balance Technology (OBT) and provides a long-lasting hydration for the skin. Results with Restylane Volyme Lidocaine are immediate visible after the first treatment. In general, results depending on the individuals lifestyle, age and technique used by the medical professional. Restylane Volyme contains lidocaine for a more comfortable injection.

CLINICAL STUDIES:

High satisfaction has been demonstrated with more than 90% wanting to have the treatment again and more than 95% stating the result looked natural. Treatments with Restylane Volyme offer a particularly high level of comfort due to the integrated local anaesthetic, lidocaine.

Restylane Volyme is safe to use and produces instant, long-lasting results using a non-surgical method of injection. As hyaluronic acid is a naturally occurring substance present in all living organisms. For the safest results, the manufacturer cautions patients to avoid exposure to extreme temperatures after treatment until the expected inflammation has disappeared.



Injectable Gel With Lidocaine

SKINBOOSTERS™ VITAL

What is Restylane SkinBoosters with Lidocaine?

Restylane Vital Lidocaine is a skinboosters designed to improve the skin's quality by boosting the hydration levels deep inside the skin and increasing its smoothness, elasticity and firmness. Restylane Vital Lidocaine is suitable for delicate skin and is indicated for areas, like the Face, lips, neck, décolletage, back of hands, skin around knee, skin folds armpits and elbow skin. Restylane Vital Lidocaine contains lidocaine for a more comfortable injection. The recommended number of treatments are three or more injections, each administered every 4 weeks.

Launched in 2004 by **Galderma** laboratories, Restylane Skinboosters is the pioneering treatment of hyaluronic acid for deep hydration.

Some 60 published clinical studies have highlighted its effectiveness as well as the sustainability of the results obtained (up to 12 months).

The smart click system a built in dose guide creates a clicking sound with each injected dose. Activate the smart click system by locking the button located on the finger grip into place.

What are the Benefits?

Restylane Skinboosters are excellent for improving skin texture and tone. Unknown to most, Skinboosters also improve acne scars by plumping them up. The benefits of the treatment include:

- Smoother skin
- Improvement in fine lines and wrinkles
- Firmer skin (face & body)
- Brighter, more luminous skin (face & body)
- Better skin hydration (face & body)
- Improvement in acne scars



Restylane SkinBoosters

Refresh Your Skin



**Restore
Your Beauty
With Restylane Lyft**

Injectable Gel With Lidocaine

LYFT™

The PERLANE Collection

What is Restylane Lyft with Lidocaine?

Restylane Lyft with Lidocaine is a clear injectable gel composed of hyaluronic acid, a natural substance that already exists in the body. Restylane Lyft with Lidocaine is non animal-based and free from animal protein. Allergy pretesting is not necessary. The lidocaine in Restylane Lyft with Lidocaine has been added to reduce the discomfort associated with the treatment.

CLINICAL STUDIES:

Restylane® Lyft works in two phases. While some results are visible immediately, the fillers do take some time to reach their full effect. After the injections, the gel adds volume that immediately improves the contour of face or hands.

A clinical study found that 95 percent of patients showed improvement within two months of treatment, and over half of the approximately 200 participants maintained the fullness for over a year or longer.

High satisfaction is demonstrated with more than 90% satisfied with how attractive balanced their face looked after a Restylane Lyft treatment.



Find your natural looking
beauty here



97%
satisfied

were satisfied
with their treatment results

With over 40 million treatments
worldwide







95%
satisfied

would recommend
the products to a friend

Restylane
GALDERMA



Product	HA Concentration	Durability	Volume	Indication	Level of Injection
	HA:20mg/ml Lidocaine: 0.3 % OBT	Up to 18 months	1ml	- shape and contours	 Hypoderm and bone contact
	HA:20mg/ml Lidocaine: 0.3 % NASHA	Up to 12 months	1ml	- shape and contours	 Hypoderm and bone contact
	HA:20mg/ml Lidocaine: 0.3 % OBT	Up to 12 months	1ml	- filling line and wrinkles deep -voluming	 Hypoderm and bone contact
	HA:20mg/ml Lidocaine: 0.3 % Non cross-linked NASHA	Up to 12 months	1ml	- face - back of hands - neck - lips - décolletage - skin around knee - skinfolds armpits - elbow skin	 Subcutaneous